

Passion to innovate An introduction from our President

Fantastic five				
Unique 3D combination – an industry first				
CBCT				
3D face photo10				
3D model scanning12				
Real-time jaw movement – in 3D14				
Planmeca ProMax [®] 3D family Key features16				
Ease of operation18				
Advanced technology20				
Pioneering low dose 3D imaging22				
2D and 3D imaging with one sensor24				
Extraoral bitewings26				
Quality cephalometry for orthodontics				
Planmeca Romexis [®] – one software for all your needs				
Reinventing 3D imaging32				
The complete implant workflow				
3D tools for orthodontists and dental labs				
Access to unique X-ray device data				
Your mobile world of imaging				
Share images and expertise online40				
Professionals proudly present the Planmeca $ProMax^{\circ}$ 3D family				
Planmeca ProMax [®] 3D s44				
Planmeca ProMax [®] 3D Classic46				
Planmeca ProMax [®] 3D Plus48				
Planmeca ProMax [®] 3D Mid50				
Planmeca ProMax [®] 3D Max52				
Patient specific implants				
Stand out with colour				
Technical specifications				



I'm extremely proud of our product innovations, and for over 40 years we've worked closely with dental professionals to set new standards in our field. What makes us a bit different is that all core product development and manufacturing takes place in Finland – ensuring exceptional quality and unmatched attention to detail at every stage of the process.

Heikki Kyöstilä President and founder Planmeca Group



"Welcome to the future of digital imaging. It gives me great pleasure to introduce you to our world-leading Planmeca ProMax[®] 3D X-ray units and Planmeca Romexis[®] imaging software – with a pioneering combination of 3D images that takes you closer for an even greater understanding of what your patients need.

This brings us to our **Planmeca ProMax**[®] product family, taking care of all your 2D and 3D imaging needs in a single unit. Each product is a true all-in-one unit, offering easy-to-use controls and incredible patient comfort. We have a dedicated team of in-house R&D professionals behind the scenes, all determined to make the best possible products for you and your patients. Therefore I am thrilled to invite you to discover our complete selection of advanced 3D solutions."



Fantastic five

Meet the Planmeca ProMax[®] 3D family



Planmeca ProMax[®] **3D** is a product family consisting of exceptional all-in-one units. With three different types of three-dimensional imaging – as well as panoramic, extraoral *bitewing and cephalometric imaging – these intelligent* products can meet all your maxillofacial imaging needs.

True all-in-one units for all your imaging needs.

Mac and Windows compatible

Unique 3D combination – an industry first

We're the first company to combine three different types of 3D data with one X-ray unit. The **Planmeca ProMax**[®] **3D** family brings together a Cone Beam Computed Tomography (CBCT) image, 3D face photo and 3D model scan into one 3D image – using the same advanced software. This 3D combination creates a virtual patient in 3D, helping you with all your clinical needs.






CBCT Why Cone Beam Computed Tomography?

Cone Beam Computed Tomography (CBCT) is an X-ray imaging technology where a large number of 2D images are taken of a patient from different angles. A 3D volumetric image is then calculated from these 2D projections. The resulting images can be viewed with our advanced imaging software from any angle, including the axial, coronal, sagittal and cross-sectional planes.

See more than ever before





Renowned dental implant surgeon Franck Renouard couldn't imagine working without his CBCT

"I acquired my Planmeca ProMax® 3D Classic in 2007 and was one of the very first users in France. The choice was quickly made, as Planmeca's unit was far ahead of its competitors."

All necessary diagnostic information from a single unit

"In implant cases, I usually start the analysis by taking a panoramic radiograph or a simple intraoral radiograph. As soon as I discover an ambiguity or low volume, I go for CBCT. For some indications, such as sinus lifts or onlay bone graftings, I always use CBCT.

I usually take large-volume studies right away. This enables diagnosis of endodontic or bone lesions that could otherwise go undetected in sites other than the implant site. The unit's resolution is more than sufficient for everyday examinations. In cases where a patient has a lot of metal restorations in their mouth or a problem with staying still, the artefact filter is very useful."

Full sinus visualisation

"I take a CBCT study systematically before a sinus lift. It provides visualisation of the sinus anatomy and allows me to see if there is an intrasinus pathology or anatomical features such as septa. It also allows detecting possible thick antral arteries, which are common in the bone wall. I need to be aware of these parameters before surgery.

Nowadays, I do not like to receive paper-based exams, as radiologists often do not provide the axial views which are essential in sinus study. When I take the radiographs myself, I can choose the slice that interests me."

Dr Renouard specifies that he has not received any financial compensation or other benefit for this







Improved patient satisfaction

I then use CBCT to check the integration of my graft when I fill the tissue before implant placement. I also need CBCT when there are complications, which occur in 3-5% of cases. I can assure you that the patients are delighted to receive their diagnosis and treatment immediately, with no need to visit a specialist radiology centre.

Planmeca ProMax 3D Classic is a well-designed radiology unit. All in all, using 3D has become natural to me. I only wonder how we managed before 3D!"

Dr Franck Renouard, DDS, Paris, France

Unique 3D combination

3D face photo

Planmeca ProFace[®] is an exclusive 3D face photo system available for all of our 3D X-ray units. This pioneering integrated system produces a realistic 3D face photo and CBCT image in a single imaging session. You can also take a separate 3D face photo without exposing your patient to any radiation.

The world's first X-ray integrated face camera



Designed to fulfil the most diverse diagnostic needs of today's maxillofacial and dental professionals, **Planmeca ProFace**[®] is a highly effective tool for pre-operative planning and treatment follow-up. It's also ideal for patient motivation and for sharing information with colleagues.


Safer and faster facial surgery

The 3D photo visualises soft tissue in relation to dentine and facial bones. As both a CBCT image and a 3D photo are generated in one imaging session, the patient position, facial expression, and muscle position remain unchanged – resulting in images that are perfectly compatible.

Careful pre-operative planning – where you can study the facial anatomy thoroughly using our **Planmeca Romexis**[®] software – facilitates accurate and detailed operations and enhances the aesthetic result.









Pre and post-operative comparisons

Measure distances and relationships between bone and soft tissue


Superimpose images for comparison


Shows deviation to instantly see any changes

Create a 2D photo

series automatically









Unique 3D combination

3D model scanning

You can use all X-ray units in the Planmeca ProMax[®] 3D family to scan both impressions and plaster casts – an exciting feature that was an industry first for our CBCT units. And with our advanced Planmeca Romexis[®] software, the digitised models are available immediately and stored for later use.



Scanning a plaster cast to a digital model



Scanning an impression to a digital model

Advantages of 3D model scanning

Digital models save space


3D digital models are stored in the Planmeca Romexis® database in standard STL format, which reduces the need to make or maintain physical plaster casts.

Create your virtual patient


The scanned 3D model can be superimposed on to CBCT data, creating a virtual patient and helping you with all your clinical and treatment planning needs. The combined data set provides an artefact-free model of your patient's dentition including bone, crowns and soft tissue. This offers valuable new options for implant planning, surgical guide manufacturing, orthodontic purposes and orthognathic surgery.












Scanned impressions of upper and lower arches and bite index in 3D

Upper and lower arch models in occlusion. A useful tool for orthodontic treatment planning and patient progress follow-up

Superimposed CBCT and 3D model of upper jaw. Measure, compare,





Crown, impression scan, and CBCT for more accurate implant planning

Analyse the STL data further in the Planmeca Romexis® 3D Ortho Studio module, and then carry out a comprehensive dental cast analysis and create an orthodontic treatment plan

Unique 3D combination

Real-time jaw movement – in 3D

Planmeca 4D[™] Jaw Motion is the only true CBCT integrated solution for tracking, recording, visualising and analysing jaw movement in 3D. It offers incomparable visualisation and measurement data of mandibular 3D movements in real-time – creating a fourth dimension in diagnostics.



Key components of Planmeca 4D[™] Jaw Motion

- CBCT image of a patient, for example a Planmeca Ultra Low Dose™ image
- Planmeca ProMax[®] 3D Mid or Planmeca ProMax[®] 3D Max X-ray unit equipped with the Planmeca ProFace[®] face photo option
- Planmeca Romexis[®] 4D Jaw Motion software module
- Special glasses and a lower jaw tracking device with light-weight reflective spheres

Applications include:

- Temporomandibular disorder (TMD) diagnostics
- Mandibular movement analysis
- Articulator programming
- Condyle-Fossa relationship during jaw movement
- Preoperative planning
- Postoperative treatment verification

Key features:

- The only CBCT integrated jaw tracking solution
- Track, visualise and record real-time jaw moveme
- Visualise movements in the Planmeca Romexis®
- Record movements for later use and analysis
- Measure and visualise the movement path of on
- Export movement and measurement informatio format for analyses and treatment planning
- Align digital dental models with a CBCT image for

ent in 3D
software without delay
e or more points of interest in a 3D image
on to 3 rd party software in an XML or CSV
or occlusion analysis



Planmeca ProMax[®] 3D family Key features

Advanced technology:

- Ideal resolutions and optimal balance between image quality and patient dose always complying with the ALARA (*As Low As Reasonably Achievable*) principle
- The pioneering Planmeca Ultra Low Dose[™] protocol enables CBCT imaging with an even lower dose than traditional 2D panoramic imaging
- Optimal volume size and location for every clinical need
- Special imaging protocols for dental and ENT applications
- Certified for use with the **suresmile** system for orthodontics

Effortless use:


- Effortless patient positioning and unmatched comfort
- True all-in-one X-ray units not only for 3D imaging, but 2D panoramic and cephalometric imaging as well
- Easy to use for a smooth workflow
- Planmeca Romexis[®] software
- Mac and Windows support

Some modalities mentioned may not be compatible with the full range of X-ray units in the Planmeca ProMax[®] 3D family. For full details on availability, see technical specifications at the end of the brochure.



Ease of operation

Our **Planmeca ProMax**[®] **3D** units are known across the world for incredible ease of use and exceptional patient comfort. A relaxed patient means a smooth imaging workflow and the best quality images.



Open patient positioning

- Effortless positioning with open-face architecture
- Unrestricted view of your patient
- No claustrophobic feeling for your patient
- Fine adjustment using positioning lasers and joystick
- Verify correct positioning with a scout image

Scout images for easy

Scout images and 2D views help

positioning and can even be used for

positioning

preliminary diagnosis.


• Easy wheelchair accomodation with side-entry access

User-friendly Planmeca ProTouch[™] control panel

- Clear and straightforward graphical user interface guides you smoothly through the work process
- Pre-programmed sites and exposure values for different image types and targets save you time and allow you to focus on your patients

Easy imaging with ready-designed protocols

- Imaging protocols designed for specific diagnostic tasks, areas, or target sizes
- Appropriate volume size, resolution, and exposure values
- Automatic selection and adjustment of the target position
- Reduced volume sizes for child patients to prevent unnecessary radiation







Advanced technology

Our intelligent high-tech solutions and algorithms guarantee an ideal imaging geometry, perfect usability, and crystal-clear images free from noise and artefacts.

> Intelligent solutions for the best image quality





Original volume, 400 µm



Reconstructed volume with higher resolution, 200 µm

SCARA technology

The precise, free-flowing, computer-controlled SCARA (Selectively Compliant Articulated Robot Arm) arm construction can produce any movement pattern required. This enables accurate and reliable volume positioning and volume diameter adjustment, reducing the amount of radiation your patients are exposed to.

New 120 kV tube voltage

120 kV tube voltage enables optimised image quality for challenging targets - reducing artefacts and ensuring higher contrast images.

Optimised imaging modes for different needs

- Low dose mode captures an image with a minimal dose of radiation. Ideally suited for orthodontic, pediatric and sinus studies. Voxel size 400 or 600 µm
- Normal mode is the best choice for most common imaging needs. Voxel size 200 µm
- High definition mode is designed for imaging of small objects, such as ear bones. Voxel size 150 µm
- Braces protocol offers optimised exposure settings for imaging patients with brackets. Voxel size 150 µm
- High resolution provides more detail, when necessary. Voxel size 100 µm
- Endodontic mode offers the best resolution with the smallest size. Voxel size 75 µm

ROI for higher resolution images

The ROI (Region of Interest) reconstruction function can generate a new small voxel volume from the image data of a previously taken large voxel volume. This enables a more precise diagnosis without the need for an additional radiation dose for the patient.


New endodontic mode

The endodontic imaging mode provides perfect visualisation of even the finest anatomical details. This advanced imaging mode is an ideal choice for endodontics and other cases with small details.

- Extremely high resolution with 75 µm voxel size
- · Enables precise diagnostics and treatment plans

Noise-free images with Planmeca AINO[™]

Planmeca AINO[™] is an intelligent 3D noise filter that removes noise from CBCT images without losing valuable details. The filter is useful in the Planmeca Ultra Low Dose™ protocol, where noise is induced by the particularly low dose. It also improves image quality in the endodontic imaging mode, where noise is inherent due to an extremely small voxel size. Planmeca AINO also allows to reduce exposure values in all other imaging modes by reducing noise.




Without noise removal

Planmeca ARA[™] – enjoy your **3D images without artefacts**


Metal restorations and root fillings in the patient's mouth can cause shadows and streaks in CBCT images. The intelligent **Planmeca ARA**[™] system removes these artefacts efficiently from Planmeca ProMax® 3D images. Our pioneering algorithm is based on a heuristic model and is the result of extensive scientific research and a vast amount of clinical patient data.



Without artefact removal







With Planmeca AINO[™] noise filter



With Planmeca ARA[™] artefact removal algotrithm

Pioneering low dose 3D i maging

Planmeca ProMax[®] 3D units offer a unique **Planmeca Ultra Low Dose[™]** imaging protocol that enables CBCT imaging with an even lower patient radiation dose than standard 2D panoramic imaging.

More information, less radiation

Planmeca Ultra Low Dose[™] can be used with all voxel sizes and in all imaging modes from Normal to Endodontic mode. Using the Planmeca Ultra Low Dose protocol reduces the effective patient dose by an average of 77% without a statistical reduction in image quality^{*}.

The unique and pioneering imaging protocol is based on intelligent 3D algorithms developed by Planmeca. Our 3D imaging system always allows the clinician to choose the optimal balance between image quality and dose, based on the ALARA principle.


* Study of Orthodontic Diagnostic FOVs Using Low Dose CBCT protocol (Ludlow, John Barrett and Koivisto, Juha). planmeca.com/ULD-poster


Ideal for many clinical cases

The Planmeca Ultra Low Dose protocol has proven to be ideal for many clinical cases.


- Orthodontics:
- Defining the amount of bone around the root
- Localising unerupted and impacted teeth before orthodontic treatment
 Defining orthodontic landmarks for
- cephalometric analysis
- Post-operative and follow-up images in maxillofacial surgery
- Airway studies
- Sinus studies
- Implant planning












Planmeca ProMax[®] 3D Classic

- + FOV Ø 40 x 50 mm / Voxel size 150 μm
- Effective patient dose 14.4 µSv



Planmeca ProMax[®] 3D Max

- + FOV Ø 85 x 50 mm / Voxel size 400 μm
- Effective patient dose 4.0 µSv

The Planmeca Ultra Low Dose[™] protoco has changed 3D imaging completely

We at MESANTIS[®] 3D DENTAL-RADIOLOGICUM produce about 7,500 CBCT images per year at eight locations in Germany.

Our main concern in X-ray imaging is to reduce the possible radiation dose as much as is reasonably achievable (ALARA principle). Traditional digital 2D X-rays at an orthodontist's clinic usually have an effective dose ranging between $26-35 \ \mu$ Sv (ICRP 2007). Conventional CBCT images of the head with modern CBCT equipment show an effective dose ranging between $49-90 \ \mu$ Sv.

The latest image protocol with a specific associated algorithm is called the Planmeca Ultra Low Dose[™] protocol. In medical terms, it allows radiologists to adjust imaging parameters individually according to the clinical needs of each case. The mA-values, in particular, can be individually adjusted and reduced for each patient, as it is required according to all international scientific guidelines. Therefore, it is possible to further reduce the effective dose significantly by using the Planmeca Ultra Low Dose protocol. Depending on the field of view, nowadays CBCT equipment with a Planmeca Ultra Low Dose algorithm has an effective dose between 4 to 22 or 10 to 36μ Sv.

Our patients and referring colleagues are always happy to hear that the effective dose for certain indications is now even lower than in traditional 2D X-ray imaging. Since last year we have been able to replace the common CBCT protocols with the Planmeca Ultra Low Dose protocol.

At MESANTIS[®] 3D DENTAL-RADIOLOGICUM in Germany, the Planmeca Ultra Low Dose imaging protocol is used either with a small or large field of view. Using the new protocol, a lot of patients can benefit from improved 3D diagnostics without being exposed to a higher radiation dose.

Prof. Dr. Axel Bumann

Prof. Dr. Axel Bumann

DDS, PhD, Orthodontist,

Oral surgeon, Oral and

Maxillofacial Radiology,

DENTAL-RADIOLOGICUM

MESANTIS[®] 3D





Planmeca ProMax[®] 3D Mid

- FOV Ø 200 x 170 mm / Voxel size 600 μm
- + Effective patient dose 14.7 μSv





Planmeca ProMax[®] 3D Mid

- FOV Ø 200 x 170 mm / Voxel size 600 μm
- Effective patient dose 29.2 μ SV

2D and 3D imaging with one sensor

Our advanced **SmartPan**[™] imaging system uses the same 3D sensor also for 2D panoramic imaging.




2D SmartPan[™] – unique panoramic imaging


- A unique system for 2D imaging
- Uses the same 3D sensor for 2D panoramic imaging, eliminating the need to change sensors
- Users can browse between panoramic images and select the most suitable one for diagnosis
- Same patient positioning and image processing parameters as in 2D imaging programs

2D programs

Standard:	Standard panoramic
Basic panoramic programs	Lateral TMJ (closed & open)
	PA TMJ (closed & open)
	PA sinus
Standard	Child (Paediatric) mode for each standard and optional program to reduce the dose
Optional	Horizontal and vertical segmenting for panoramic program
Optional	True Bitewing
<i>Optional:</i> Advanced	Interproximal panoramic
panoramic programs	Orthogonal (perio) panoramic
	Lateral-PA TMJ
	Lateral multiangle TMJ
	PA multiangle TMJ
	PA linear sinus
	Lateral sinus
L	1



Normal **SmartPan**[™] produces 9 different parallel panoramic layers with about 2 mm shift and one autofocus layer



MultiView SmartPan[™] calculates 9 different rotated panoramic layers. This allows adjusting the view angle for improved diagnosis.



Extraoral bitewings

What if you could do all your routine diagnostic imaging extraorally?


Planmeca ProMax[®] extraoral bitewings are ideal for periodontics, elderly and child patients, claustrophobic patients, patients with a strong gag reflex, and patients in pain. Extraoral bitewings enhance clinical efficiency and take less time and effort than conventional intraoral bitewing imaging.

What are the advantages of extraoral bitewings?

- Ideal for all patients no sens positioning required
- Consistently opens interproximal contacts, giving better diagnostic value
- Larger diagnostic area than in intraoral modalities
- More clinical data: canine to third molar
- Enhanced clinical efficiency takes less time and effort than conventional intraoral bitewings
- Enhanced patient experience and comfort eliminates gagging

Better diagnostic value with extraoral bitewings








Standard panoramic image of the same patient as the bitewing above





True Bitewing program, 5-year-old child



True Bitewing program, 8-year-old child



True bitewings only possible with our SCARA3 technology

Quality cephalometry for orthodontics

We offer exceptional equipment and the most advanced software for all your orthodontic needs.

Cephalometric imaging with Planmeca ProMax[®] units

- The functional and easy-to-use head positioner ensures accurate positioning for all cephalometric projections
- The carbon fibre ear posts and nasal positioner are extremely stable, hygienic, and transparent to radiation
- The unit automatically aligns itself to take cephalometric exposures and then selects a corresponding collimator
- The rotating tube head in the 3D unit eliminates the need to remove the 3D sensor

Two available options:

One-shot Planmeca ProCeph™ cephalostat

- · Effective one-shot cephalostat
- · Short exposure time no motion artefacts, low patient dose
- Image sizes from 18 x 25 cm to 30 x 25 cm

Easier and more accurate than ever before

350

Scanning Planmeca ProMax[®] cephalostat

- Digital cephalostat that scans your patient's head horizontally using a narrow X-ray beam with an extremely low effective dose of radiation
- · Exceptional flexibility in image formats, with field sizes of up to 30 x 27 cm



Planmeca Romexis[®] **Cephalometric Analysis module**

- Create cephalometric analyses and superimpositions in minutes
- Fully customisable analyses, norms and reports
- Microsoft Excel export and import function
- Compatible with Windows operating system

Planmeca Romexis® one software for all your needs

We offer a revolutionary all-in-one software solution for clinics of all sizes. Our world-leading **Planmeca Romexis**[®] software is the brains behind all of our products, bringing together all the devices at your dental clinic from CAD/CAM to imaging devices and dental units. It supports the most versatile range of 2D and 3D imaging modalities.

Planmeca Romexis

2D imaging

Infection contro

information and monitoring

Dental units


Imaging and CAD/CAM in one software – an industry first



Planmeca Romexis®

Reinventing 3D imaging

Our pioneering **Planmeca Romexis**[®] software offers specially designed tools for implantologists, endodontists, periodontists, prosthodontists, orthodontists, maxillofacial surgeons, and radiologists. You can also view your images wherever you are using our mobile apps, and enjoy unmatched compatibility with other systems.



Excellent tools for quality images

With a complete set of tools for image viewing, enhancement, measurement, drawing and annotations, **Planmeca Romexis**[®] improves the diagnostic value of radiographs. Versatile printing and image import and export functionalities are also included. The software consists of different modules – so you can choose those most suited to your needs.


Convenient 3D diagnosis

The Planmeca Romexis 3D rendering view gives an immediate overview of the anatomy and serves as an excellent patient education tool. The images can be instantly viewed from different projections or converted into panoramic images and cross-sectional slices. Measuring and annotation tools – such as nerve canal tracing – assist in safe and accurate treatment planning. Free Planmeca Romexis[®] Viewer application planmeca.com/Viewer

> Full-featured viewer application No installation required Mac and Windows support Distribute to specialists or patients

Superimpose CBCT

New to Planmeca Romexis 3D, the module allows the superimposition of two CBCT images. It is a valuable tool for before-and-after comparisons and can be used for orthognathic surgery follow-ups, as well as orthodontic treatments, for example. The module also allows users to compare CBCT and MRI images side by side – providing a comprehensive view of a patient's anatomy.





Tooth segmentation

Planmeca Romexis provides a new, intuitive and efficient tool for segmenting a tooth and its root from a CBCT image. The guided process enables quick segmentation of a patient's full dentition. Surface models of segmented teeth can be visualised, measured and utilised e.g. in **Planmeca Romexis® 3D Ortho Studio** as part of orthodontic treatments.


Easy sharing of results

Studies can be quickly converted into multi-page printouts or handed out with the free **Planmeca Romexis® Viewer** media. Cases can be seamlessly transferred to mobile devices or partner clinics that also use Planmeca Romexis.

Best compatibility with other systems

Planmeca Romexis offers excellent compatibility with other systems, allowing you to freely use third-party products at your clinic. TWAIN support and DICOM standard compliance ensure that our flexible software can be used effortlessly with most systems.





Visualise and measure airways and sinus volumes before and after treatment for simplified diagnosis and treatment planning. Our advanced software tools allow accurate measurements in 3D space. Measurements can easily be reviewed using the saved views. Planmeca Romexis[®]


The complete implant workflow

Our **Planmeca Romexis**[®] **3D Implant Planning** module offers all the necessary tools for fully digital implantology – from planning to guided surgery. The software's implant library includes realistic implant models and abutments, as well as collections of sleeves for guided surgery. After completing the implant plan, a surgical guide can be immediately designed in the same **Planmeca Romexis**[®] software with just a few clicks.



Realistic implant models from over 60 manufacturers

The **Planmeca Romexis**[®] software platform provides the perfect environment for top-down implant planning. By superimposing a crown and dental model onto CBCT data, users can create a complete virtual setup for optimally positioning the implant – taking prosthodontic and surgical perspectives into account.





Mark the nerve on the CBCT image

Superimpose the 3D model scan onto the CBCT image with the Planmeca Romexis[®] software

Use the Planmeca Romexis[®] crown library, or import a patient-specific crown from the CAD system to the software


Select the preferred implant (abutment) and sleeve from the extensive Planmeca Romexis[®] library and find the optimal position for it from a prosthetic and surgical perspective

Design the surgical implant guide with just a few clicks in Planmeca Romexis[®] – the software will create an open STL file of the design


Print the surgical guide with Planmeca Creo™ or any other suitable 3D printer Planmeca Romexis®


3D tools for orthodontists and dental labs

Planmeca Romexis[®] 3D Ortho Studio provides innovative tools for orthodontists and dental laboratories. Our advanced module is designed for the examination and analysis of digital dental models scanned with Planmeca ProMax[®] 3D X-ray units or the Planmeca PlanScan[®] intraoral scanner – and also for planning orthodontic treatments in 3D.















Dental model analysis

Dental impressions and plaster casts scanned with the Planmeca ProMax[®] 3D model scanning mode can be aligned with the bite index using the Planmeca Romexis[®] software. Examination, analysis and treatment planning are then conveniently done in the Planmeca Romexis® 3D Ortho Studio module.

The module makes dental model analysis easier then ever by offering all the necessary tools for virtual base creation, occlusion inspection, and versatile teeth and arch measurements.





Planmeca PlanScan



Plaster cast in Planmeca ProMax[®] 3D




Impression scan in Planmeca ProMax® 3D

Treatment planning in 3D

A staged treatment plan can be established in Planmeca Romexis 3D Ortho Studio by displacing the teeth in a virtual tooth setup, while visualising intersections and contacts.

All the applied changes such as tooth movements, interproximal reductions and tooth extractions are summarised in a detailed treatment plan report. The treatment plan can be easily shared with others.



Export of digital dental models in STL format

Planmeca Romexis 3D Ortho Studio generates a sequence of digital dental models for each treatment stage. The models can be exported in STL format for 3D printing and custom appliance design and manufacturing.

The module is compatible with the Windows operating system









The Braces imaging protocol for Planmeca ProMax[®] 3D units

Our special Braces* imaging protocol offers optimised exposure settings for imaging patients with brackets. The acquired images can be utilised when designing braces. The Braces imaging protocol is optimised for use with suresmile.

suresmile

to be sure.

Suresmile[®] certification

Planmeca ProMax 3D units have been certified for use with suresmile by OraMetrix.

The comprehensive treatment management system suresmile is designed to enable orthodontists to visualise and simulate multiple diagnostic set-ups and design customised archwires for every patient. The accuracy of the patient scans plays a critical role in maximising the effectiveness of the system.

*Available for Planmeca ProMax® 3D Classic, Planmeca ProMax® 3D Mid and Planmeca ProMax® 3D Max.

Access to unique X-ray device data

Our **Planmeca Romexis[®] Clinic Management** software module provides valuable and reliable equipment usage data. It is an ideal solution for all – from smaller private practices to large multi-site clinics.

All Planmeca CBCT and panoramic units can be connected to the Planmeca Romexis[®] Clinic Management software module. The module collects information on the number and type of taken images and provides data from different Planmeca ProMax[®] and Planmeca ProOne[®] X-ray units for reporting, improving clinic efficiency, and planning new investments.

Planmeca Romexis Clinic Management helps users stay on top of clinic operations with real-time information on Planmeca CBCT and panoramic X-ray units and their status. Detailed equipment information – such as software versions, calibration dates, and history logs - ensures high-quality service and patient care at all times.

Your mobile world ofimaging Planmeca mRomexis[™] image viewing application

Our advanced **Planmeca mRomexis**[™] multiplatform application allows you to flexibly access your images on the go. Remove the constraint of place – consult with colleagues and communicate with patients easily wherever you are.

Stay mobile with the Planmeca mRomexis[™] image viewing application

Our fast, easy and light **Planmeca mRomexis**[™] image viewing application is designed for flexible multiplatform use. It is available for iOS and Android devices and as a browser-based desktop application. Access all your images in the Planmeca Romexis® database on a local network or carry images with you on your mobile device. Experience the new level of freedom Planmeca's mobile world can offer!





NETWORI

26.5

main 27 and





Planmeca mRomexis allows you to stay informed at all times. It is the ideal solution for fast and light viewing of 2D and 3D images, 3D surface models and Planmeca ProFace® facial photos. Images can be flexibly shared via the Planmeca Romexis[®] Cloud image transfer service.

Constantly keep up with your workflow with Planmeca mRomexis and ensure that images most relevant to your next task are always readily available.





Download the Planmeca mRomexis[™] application for iOS and Android from the App Store or Google Play.



Planmeca Romexis[®]

Share images and expertise online

Planmeca Romexis[®] Cloud is a secure image transfer service for **Planmeca Romexis[®]** users and their partners. Now you can easily share images and CAD/CAM cases with any specialist or patient.

Planmeca Romexis[®] user

- Radiology center
- General practice

Planmeca Romexis[®] Cloud IMAGES REFERRALS INTERPRETATIONS TREATMENT PLANS

Advantages

- Seamlessly integrated into
 Planmeca Romexis[®] ensuring an efficient workflow – no need for external applications or CDs and DVDs
- Automatic delivery of images and attachments
- Automatic notification to recipient of new cases

Features

• Cases can be sent to any recipient who

• Streamline your communication with

has an e-mail account

information

• Secure transfer and storage of

Planmeca Romexis® Cloud

- Sending images to recipient
- 2D images: panoramic, cephalometric, photos, intraoral X-ray images
- 3D images: CBCT, 3D photos, surface scans
- All annotations and other elements are included

Sending documents to recipient

 Attach one or more referrals, reports, or other documents

Versatile possibilites for communication

Recipients can download and view images at no cost using:

- Planmeca Romexis
- Planmeca mRomexis[™] image viewing application for iOS and Android

Free Planmeca Romexis[®] Viewer

Planmeca Romexis[®] software and

Planmeca Romexis® Cloud subscription are required for sending new cases. Visit http://online.planmeca.com/ to subscribe and start sending images now.





Anybody, anywhere

- General practitioner
- Colleague
- Radiologist
- Specialist
- Dental lab
- Patient

Professionals proudly present the Planmeca ProMax[®] 3D family



Which one is rig

Planmeca ProMax[®] 3D s

Planmeca ProMax[®] 3D s is an ideal 3D un It is perfect for single implant, endodont

Planmeca ProMax[®] 3D Classic

The Planmeca ProMax[®] 3D Classic imagin area, so the unit gives a clear view of the

Planmeca ProMax[®] 3D Plus

The newest member in our 3D family, Pla variety of different volume sizes and is a

Planmeca ProMax[®] 3D Mid

Thanks to its wide volume size selection, a wide range of diagnostic tasks without

Planmeca ProMax[®] 3D Max


Planmeca ProMax[®] 3D Max is a dedicated all required volume sizes when diagnosir - from the smallest special cases to imag


The interviewed have not received any financial compensation or other benefit for the interviews that follow.

ht for you?	• • • •
it for capturing small details. ic, and wisdom tooth cases.	• • • •
ng sensor covers the whole dentition mandible and maxilla.	• • • •
nmeca ProMax [®] 3D Plus, offers a wide great choice for any imaging need.	• • • •
Planmeca ProMax [®] 3D Mid handles compromising best practices.	• • • •
d 3D imaging device that produces ng the maxillofacial region res of the entire head.	• • • •



Planmeca ProMax[®] 3D s

















Long-term cooperation with Planmeca

"We purchased a Planmeca ProMax[®] 3D s for our dental clinic about four years ago or so. Before that, we had equipped our clinic with five Planmeca dental units, so it was only natural to continue the cooperation with Planmeca also on the X-ray side. Also, several radiologists recommended Planmeca's 3D units to us for their high quality.

We use the unit for implant cases, for lower third molar surgery, and for endodontic cases – particularly in difficult infection cases of teeth with multiple roots. Personally, I use the Planmeca Romexis® 3D Implant Planning module the most. It's very practical as I can virtually place the implants myself in the software.

The unit itself is very easy to use - our whole staff uses it, although mainly dentists take 3D images. Positioning is effortless and images are of high quality. And the unit's design is stylish and refined.

I would definitely recommend the unit to others. We have just taken the new sensor into use and I am very satisfied with the image quality. And the feedback from consulting radiologists has been good as well."

Ari Mäkelä, Licentiate in Dentistry Dental Care Center Janne, Järvenpää, Finland



Chinese hospital chose Planmeca ProMax[®] 3D s

"I bought the Planmeca ProMax" 3D s system in September 2010, so I have been using it for over 2 years now. Factors influencing my decision were Planmeca's good reputation and quality-price ratio. For me, it is also important that everyday performance is excellent and if necessary, the after sales service works quickly.


I use my Planmeca 3D s system for various cases – for diagnosis in oral and maxillofacial surgery, for implantology, for diagnosis of periodontal and dental pulp diseases, and for orthodontics. The image quality is very clear, which makes diagnosis very easy with the excellent Planmeca Romexis® software.


In implant cases, Planmeca ProMax 3D s is very important for my preparation phase. The data I get from the image of the bone structure and thickness makes the operation easy and safe for the customer.


Planmeca ProMax 3D s really adds value to my work as I can perform many different kinds of tasks quickly and efficiently."

Sun Zhizong, Dean Donggang City Stomatology Hospital, Liaoning, China

Planmeca ProMax[®] 3D Classic
















West Vantaa Dental Clinic, Finland

Finnish dental clinic chooses Planmeca ProMax[®] 3D Classic

"We decided to purchase a **Planmeca ProMax[®] 3D Classic** 8x8 for our clinic as we wanted to start taking our own CBCT images and not have to send our patients elsewhere to have their 3D X-rays taken. In such cases, there is always the risk that the treatment process will suffer due the patient's own lack of activity. Now we have our own radiologist and things have gone very smoothly. We also have two surgeons working with us, as we do a lot of implant treatments and treat also difficult endodontic cases."

Implant case acceptance has skyrocketed

"After acquiring the Planmeca ProMax 3D Classic, the amount of implant cases treated at our clinic has increased considerably. Patients are always amazed

when we offer to take their 3D images straight away. The unit is also especially suited to complicated endodontic cases, as you can notice everything in a 3D volume. It is also excellent for cases of wisdom teeth that have grown at a cumbersome angle.

The image quality produced by Planmeca ProMax 3D Classic is excellent. I think it is safe to say that we have the best 3D unit in Finland. This opinion is shared by our surgeons and many radiologists.

The **Planmeca Romexis**[®] software is a great working tool. It is logical, easy to use, and functions well – just a really good piece of software."

Pekka Nissinen, GPD, West Vantaa Dental Clinic, Finland






Optimal image quality for every single field of dentistry


"I've been using Planmeca ProMax 3D Classic ever since its introduction to the market in 2007, and have used it for all imaging purposes. The image quality has proven to be reliable in every single field of dentistry, even in the most demanding imaging cases. The unit is very user-friendly, and all in all the imaging process can be carried out in an uncomplicated manner.


The Planmeca Romexis software is, in my opinion, the best software on the market when it comes to 3D imaging."

Kim Lemberg, radiologist, West Vantaa Dental Clinic. Finland

Planmeca ProMax[®] 3D Plus













Dr. Dirk Ladig

Oral surgery practice, Hoyerswerda, Germany

German oral surgery practice is impressed with the image quality of Planmeca ProMax^o 3D Plus

"I have been using the **Planmeca** ProMax[®] 3D Plus unit in my oral surgery practice since 2013. Before that, I had good experience with Planmeca X-ray units. My panoramic X-ray unit ran smoothly for 19 years, the service was good and I was satisfied. Moreover, in 2000, I integrated cone beam computed tomography into my practice by adding a second unit. The decisive factor in purchasing the Planmeca ProMax 3D Plus unit was the radiographs of the new flat-panel devices shown to me by colleagues. The higher resolution of the images was very impressive! There was also a change in the physical layout of my practice. Instead of having two X-ray rooms, I wanted to have one. Planmeca ProMax 3D Plus combines two devices in one: OPG and CBCT. As a result, we need considerably less space.

More information in a single image

I use the device for different kinds of treatment planning; mainly implant cases, but also high-risk wisdom tooth surgery. In my view, a key benefit of the Planmeca ProMax 3D Plus is the possibility of displaying the entire mandible – including the ascending mandibular ramus and mandibular joint in a single image. I also use the images for diagnosis of foreign body location, apical variances and inflammatory processes in the jaw area. CBCT provides much better diagnostic options for screening for infectious foci in patients with unclear symptoms or certain systemic diseases. Questions related to orthodontic treatments of impacted and displaced teeth, for example, can be easily solved on behalf of colleagues.




Low radiation exposure with adjustable volume sizes


What I really like about the unit is that I can select the volume according to the required image. The radiation exposure for patients is thus kept as low as possible. I use low-dose scans particularly with orthodontic diagnosis. The layer lights are especially useful when centring the image volume.


Operating and adjusting the unit is easy. What's more, the transition from analogue to digital control went well. Since the patients stand upright within the unit, positioning them is much easier than with the predecessor of the CBCT model (with patient bench), without having any problems with motion blur. The new device is also much more pleasant for the patients because there is no feeling of constriction."

Dr. Dirk Ladig Oral surgery practice, Hoyerswerda, Germany

Planmeca ProMax[®] 3D Mid













Dr Carlo Pizzo, DDS & Dr Gioia Amico, DDS

A&P Clinic, Cittadella, Italy

Italian A&P Clinic opts for Planmeca ProMax[®] 3D Mid after a thorough market analysis

"In our new dental clinic, we have been using **Planmeca ProMax**" **3D Mid** for six months now – and we are really satisfied with it.

We chose the unit after a thorough analysis of what the market was offering. We needed an imaging unit that could provide a wide range of FOV choices, the possibility to take panoramic images and cephalometric shots, and last but not least, software that could run natively on Mac OS, because our IT infrastructure was entirely built on Apple computers. The only unit that fulfilled all of these requirements was Planmeca ProMax 3D Mid."

For every clinical application

"We love using it for taking panoramic images, preliminary treatment planning, 3D scans, wisdom teeth extractions and implant surgery. With **Planmeca Romexis**[®] – its dedicated software – we can virtually place the exact dental implants we are going to use by choosing them from the integrated 3D implant library. This feature works amazingly well."

3D magic with the latest technology


"The machine and the software work seamlessly together: they are fast, reliable and easy to use. The 3D rendering is an incredibly powerful tool for us – for visualising the real bone morphology of the patients, and for the patients themselves to understand their clinical situation and the treatment we are offering them. So Planmeca Romexis can





become a really effective communication tool. For this reason, we adopted also the **Planmeca ProFace**[®] option. By superimposing a 3D scan of the patient's face and a CBCT X-ray image, we can show our clients an easy-to-understand image, in which they can really recognize themselves. Even today, this looks like magic for many of our patients!"

Dr Carlo Pizzo & Dr Gioia Amico, A&P Clinic, Cittadella, Italy

Planmeca ProMax[®] 3D Max









Dr Corrado Gazzerro

MD, Specialist in Radiodiagnostics, **Qualified Expert in Radioprotection**

Studio Gazzerro, Genoa, Italy

Radiologist praises the versatility of Planmeca ProMax[®] 3D Max

"I was the first **Planmeca ProMax**" **3D Max** user in Italy and have been using it for about three years now. Before that, I used Planmeca ProMax[®] 3D Classic 8x8 for 2 years. And I've been using Planmeca equipment since 1995 because of their image quality, their reliability, and the fast maintenance service.

I really enjoy working with Planmeca ProMax 3D Max. I have used it for every possible dental case, including all aspects of implantology, as well as endodontics, examining alterations of the bone structure, wisdom tooth extractions,



supernumerary teeth and more. In ENT cases, I have used the unit for the study of the paranasal sinuses and facial bone structures.

One of the most remarkable advantages is the possibility to choose the image quality and therefore to optimise the patient dose. The volume selection is complete, the imaging programs are easy to use and patient positioning is effortless."

Dr Gazzerro, Studio Gazzerro, Genoa, Italy

Patient specific implants

You imagine it. We make it.

Ż

0

E.

•


Ω

Planmeca ProModel[™] offers patient-specific implants and physical models for maxillofacial surgery – all individually designed for best possible results.

> The implants are designed and manufactured to match any form, ensuring an exact fit to the patient's anatomy. The service also includes physical 3D skull models and surgical guides for assisting in both pre-planning and the surgery itself.

> > PLANMEC.

ProModel



The 3D design is created in an online meeting between surgeon and designer.

Skull model for surgical pre-planning and ready-to-use patient specific implant.

Planmeca ProModel[™] service concept

- A unique service for creating patient specific implants, surgical guides and skull models from CBCT/CT images
- 3D implants are designed in an online session between the surgeon and Planmeca designer
- Ordering is quick and easy from order to delivery in just a few business days
- Significantly lowers costs and reduces operation times by up to 4 hours
- Faster and more precise operations leading to better aesthetic results

Faster operations, precise fit and better aesthetic results





Stand out with colour

Complement the splendid design of your **Planmeca ProMax® 3D** X-ray unit by giving it a personal touch with your favourite colours. Select the perfectly matching shades from our exquisite and inspiring collection and create the looks of your dreams!



Technical specifications

Technical data

	3D s	3D Classic	3D Plus	3D Mid	3D Max
Anode voltage	60-90 kV	60-90 kV	60-90 kV	60-90 kV	60-96 kV*
				60–120 kV	60–120 kV**
Anode current	1–14 mA	1–14 mA	1–14 mA	1–14 mA	1–12 mA
Focal spot	0.5 mm, fixed anode	0.5 mm, fixed anode	0.5 mm, fixed anode	0.5 mm, fixed anode	*0.6 mm, fixed anode
					**0.5 mm, fixed anode
Image detector	Flat panel	Flat panel	Flat panel	Flat panel	Flat panel
Image acquisition	Single 200 degree rotation	Single 200 degree rotation	200 / 360 degree rotation	200 / 360 degree rotation	210 / 360 degree rotation
Scan time	7.5–27 s	9–37 s	9–33 s	9–33 s	9-40
Typical reconstruction time	2–25 s	2–25 s	2-30 s	2–55 s	2-55 s

Comparison

	3D s	3D Classic	3D Plus	3D Mid	3D Max
3D dental programs	Yes	Yes	Yes	Yes	Yes
3D ENT programs	-	-	Yes	Yes	Yes
3D face photo	Yes	Yes	Yes	Yes	Yes
3D models scan	Yes	Yes	Yes	Yes	Yes
Suresmile certification	-	Yes	-	Yes	Yes
4D jaw motion	-	-	-	Yes	Yes
2D panoramic imaging	Yes	Yes	Yes	Yes	Yes
2D cephalometric imaging	Yes	Yes	Yes	Yes	-

Dental programs

Volume size (child mode) [mm]

	3D s	3D Classic	3D Plus	3D Mid	3D Max	Voxel size, isotropic
Tooth	Ø50 x 50 (Ø42 x 42)	Ø50 x 50 (Ø42 x 42)	Ø40 x 50 (Ø34 x 42)	Ø40 x 50 (Ø34 x 42)	Ø50 x 55 (Ø42 x 50)	75 μm*, 100 μm, 150 μm, 200 μm, 400 μm
	Ø50 x 80 (Ø42 x 68)	Ø50 x 80 (Ø42 x 68)	Ø40 x 70 (Ø34 x 60)	Ø40 x 80 (Ø34 x 68)		150 µm, 200 µm, 400 µm
Teeth		Ø80 x 80 (Ø68 x 68) Ø80 x 50 (Ø68 x 42)	. ,	Ø80 x 50 (Ø68 x 42) Ø80 x 80 (Ø68 x 68) Ø100 x 60 (Ø85 x 50) Ø100 x 100 (Ø85 x 85)	Ø100 x 55 (Ø85 x 50) Ø100 x 90 (Ø85 x 75)	150 μm, 200 μm, 400 μm
	triple scan: 90 x 60 x 80	triple scan: 140 x 105 x 80				200 µm, 400 µm
Jaw			Ø160 x 50 (Ø160 x 50) Ø160 x 90 (Ø160 x 90)	Ø160 x 60 Ø160 x 100 Ø200 x 60 (Ø200 x 60) Ø200 x 100 (Ø200 x 100)	Ø230 x 60 Ø230 x 100 Ø130 x 55 (Ø110 x 50) Ø130 x 90 (Ø110 x 75)	200 μm, 400 μm, 600 μm
Face				Ø200 x 170 (Ø200 x 170)	Ø100 x 130 (Ø85 x 110) Ø130 x 130 (Ø110 x 110) Ø130 x 160 (Ø110 x 136)	200 μm, 400 μm
Skull					Ø230 x 160 Ø230 x 260	400 μm, 600 μm

ENT (Ear, Nose, Throat) programs


Volume size (child mode) [mm]


	3D Plus	3D Mid	3D Max	Voxel size, isotropic
Nose	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)	Ø100 x 90 (Ø85 x 75)	200 µm, 400 µm
Sinus	Ø90 x 90 (Ø90 x 90)	Ø100 x 100	Ø100 x 90	200 µm, 400 µm, 600 µm
	Ø160 x 90 (Ø160 x 90)	Ø100 x 170	Ø100 x 130	
		Ø160 x 100	Ø130 x 100	
		Ø160 x 170	Ø130 x 130	
		Ø200 x 100	Ø130 x 160	
		Ø200 x 170		
Middle ear	Ø40 x 50 (Ø34 x 42)	Ø40 x 50 (Ø34 x 42)	Ø50 x 55 (Ø42 x 50)	75 μm*, 100 μm, 150 μm, 200 μm , 400 μm
	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)		150 μm, 200 μm, 400 μm
Temporal bone	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)	Ø100 x 90 (Ø85 x 75)	150 μm, 200 μm
Vertebrae	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)	Ø100 x 90 (Ø85 x 75)	200 µm, 400 µm
			Ø100 x 130 (Ø85 x 110)	
Airways	Ø70 x 70 (Ø60 x 60)	Ø80 x 80 (Ø68 x 68)	Ø100 x 90 (Ø85 x 75)	200 µm, 400 µm
			Ø100 x 130 (Ø85 x 110)	
			Ø130 x 130 (Ø110 x 110)	
			Ø130 x 160 (Ø110 x 136)	

*Requires Endodontic imaging licence



Technical specifications





Dimensions

	3D s or 3D Classic	3D Plus or 3D Mid	3D Max
Α	1298-2123 mm (51.1-83.5 in.)	1315–2095 mm (51.8–82.5 in.)	-
В	1560-2385 mm (61.4-93.8 in.)	1610-2390 mm (63.4-94.1 in.)	1582-2482 mm (62.3-97.7 in.)
С	1145 mm (45.1 in.)	1130 mm (44.6 in.)	-
D	850 mm (33.5 in.)	930 mm (36.6 in.)	930 mm (36.6 in.)
E	270 mm (10.6 in.)	247 mm (9.7 in.)	222 mm (8.7 in.)
F	698 mm (27.5 in.)	810 mm (32 in.)	788 mm (31 in.)
G	1250 mm (49.2 in.)	1366 mm (53.8 in.)	1351 mm (53.2 in.)
н	777 mm (30.6 in.)	756 mm (29.8 in.)	-
J	Ø820 mm (32.3 in.)	Ø1010 mm (39.8 in.)	Ø1010 mm (39.8 in.)

Physical space requirements


	3D s or 3D Classic	3D s or 3D Classic with cephalostat	3D Plus or 3D Mid	3D Plus or 3D Mid with cephalostat	3D Max
Width	115 cm (44 in.)	200 cm (79 in.)	118 cm (47 in.)	206 cm (82 in.)	116 cm (45.3 in.)
Depth	125 cm (49 in.)	125 cm (49 in.)	137 cm (54 in.)	137 cm (54 in.)	137 cm (54 in.)
Height*	153–243 cm (60–96 in.)	153–243 cm (60–96 in.)	161–239 cm (64–94 in.)	161–239 cm (64–94 in.)	161–239 cm (64–94 in.)
Weight	113 kg (lbs 248)	128 kg (lbs 282)	131 kg (lbs 289)	146 kg (lbs 322)	131 kg (lbs 289)

Minimum operational space requirements

	3D s or 3D Classic	3D s or 3D Classic with cephalostat	3D Plus or 3D Mid	3D Plus or 3D Mid with cephalostat	3D Max
Width	150 cm (59 in.)	215 cm (85 in.)	158 cm (63 in.)	225 cm (89 in.)	158 cm (63 in.)
Depth	163 cm (64 in.)	163 cm (64 in.)	175 cm (69 in.)	175 cm (69 in.)	175 cm (69 in.)
Height*	243 cm (96 in.)	243 cm (96 in.)	239 cm (94 in.)	239 cm (94 in.)	239 cm (94 in.)

*The maximum height of the unit can be adjusted for offices with limited ceiling space.

Example installation





www.facebook.com/PlanmecaOy

Planmeca Romexis[®] imaging software

	kis intaging software
Supported 2D modalities	Intraoral
	Panoramic
	Cephalometric
	2D linear tomography
	Photos
	Stack images (CBCT slices and panoramic slices)
Supported 3D modalities	3D CBCT
SD modalities	3D photo
	3D surface scan
Supported photo sources	Intraoral camera
photo sources	Digital camera or scanner (import or TWAIN capture)
Operating systems	Win 7 Pro / Win 8.1 Pro / Win 10 (64 bit)
	Win 2008 Server / Win 2012 Server
	Mac* (OS X or newer)
	For detailed information please see system requirements of Planmeca Romexis www.planmeca.com
	*Cephalometric Analysis module, 3D Ortho Studio module and Planmeca PlanCAD Easy are supported on Windows operating systems.
Image formats	JPEG or TIFF (2D image)
	DICOM (2D and 3D image)
	STL (3D image)
	TIFF, JPEG, PNG, BMP (import/export)
Image size	2D X-ray image: 1–9 MB
	3D X-ray image: typically 50 MB–1 GB
Installation options	Client–Server
	Java Web Start deployment
DICOM 3.0 support	DICOM Import/Export
	DICOM DIR Media Storage
	DICOM Print SCU
	DICOM Storage SCU
	DICOM Worklist SCU
	DICOM Query/Retrieve
	DICOM Storage Commitment
	DICOM MPPS
Interfaces	TWAIN Client
	PMBridge (patient information and images)
	VDDS (patient information and images)
	InfoCarrier (patient information)
	Datagate (patient and user information)
3 rd party software integrations	Dolphin Imaging
	Nobel Clinician
	Materialise Dental Simplant
	Straumann coDiagnostiX
	Cybermed N-Liten







One software for all.



Planmeca Oy designs and manufactures a full line of industry-leading dental equipment, including 3D and 2D imaging devices, CAD/CAM solutions, dental care units and software. Planmeca Oy, the parent company of the Finnish Planmeca Group, is strongly committed to better care through innovation, and it is the largest privately held company in the field.

PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | fax +358 20 7795 555 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, DentroVac, Digital perfection, Economat Plus, Elegant, Flexy, Mini-dent, Perio Fresh, PlanEasyMill, Planmeca 4D, Planmeca AINO, Planmeca ARA, Planmeca CAD/CAM, Planmeca Chair, Planmeca Clarify, Planmeca Compact, Planmeca FIT, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Lumio, Planmeca Maximity, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca Romexis, Planmeca Noma, Planmeca Online, Planmeca PlanCAD, Planmeca PlanCAM, Planmeca PlanClear, Planmeca PlanMill, Planmeca Planosi, Planmeca PlanPure, Planmeca Plance, Planmeca Noma, Planmeca ProFace, Planmeca PlanCAD, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSenor, Planmeca ProX, Planmeca Xtremity, Proline Dental Stool, ProTouch, Saddle Stool, SmartTouch, Trendy and Ultra Relax are registered or non-registered trademarks of Planmeca in various countries.

