


## Daily Measure TDS (Total Dissolved Solids) of water with Water Test Pen

- DAC units water less than 3us/cm (2 ppm)
- Statims water more than 0.5us/cm (0.3ppm) but less than 10us/cm (6.5PPM)
- Prestige Autoclaves water less than 35us/cm (22ppm)
- Mammoth Autoclaves waster less than 15us/cm (9.5ppm)
- Eschamnn Autoclaves water less than 15us/cm (9.5ppm)

# **RO System Filter Schematic**





### **Carbon & Sediment Filters**

• Replace the first two vertical filters i.e sediment and carbon every 6 months.

### DI Filter

• The DI resin should be changed as the TDS rises (Assuming Carbon & Sediment changed) however it should be noted that if the TDS rises suddenly the RO membrane should be checked to ensure they are performing properly.

### RO Membrane

- It should be flushed as often as possible ideally **once per day** and at least once per week.
- The flushing of the membrane helps wash away the contaminants and salts from the membranes surface thus cleaning it, improving the membranes performance and extending it life span.
- Regular flushing extends the life of the resin.

