



Cavitron® Integrated Systems offer an open environment with the same benefits of traditional Cavitron® units.

 Compatible with Steri-Mate® 360°, a rotating handpiece that allows the clinician to fully rotate an ultrasonic insert without stopping for adjustments

 Detachable and autoclavable handpiece reduces cross contamination for enhanced infection control

 Patented Sustained Performance System, SPS™, Technology automatically adjusts power to maintain scaling efficiency, even at low power

• **BlueZone™** extended low power range designed for improved patient comfort





Cavitron Integrated Systems integrate easily with A-dec® 500 and A-dec® 300 delivery systems.

For dental professionals searching for the streamlined look of integration, you are now able to combine Cavitron Systems with the convenience of integration into the A-dec 500 or A-dec 300 dental units. You can also put more technology at your fingertips with the A-dec deluxe touch pad, which gives you an intuitive display, easy-touch control, and programmable memory for preferred ultrasonic power settings.

Item	A-dec Integrated Part #	Description	Compatible With:
Cavitron A-dec Built-In Ultrasonic Scaler Kit	54.0599.00	Cavitron Steri-Mate 360° Kit, 61"	A-dec 532, 533, 542, 332, 333, 334, 335, 336
	54.0597.00	Cavitron Steri-Mate Kit, 61"	
	54.0600.00	Cavitron Steri-Mate 360° Kit, 82"	A-dec 541
	54.0598.00	Cavitron Steri-Mate Kit, 82"	

For an integrated Cavitron kit contact:

1.800.547.1883 Within USA 1.503.538.7478 Outside USA/Canada www. a-dec.com



For a built-in Cavitron kit for field installation contact:

DENTSPLY Professional 1.800.989.8826 www.dentsply.com

