



MULTIPLE SURGERY ASPIRATORS

# TURBO SMART 2V

## Control More Power

Cattani 'Turbo Smart 2V' is the first direct replacement for any wet-ring system, suitable for up to 4 surgeries producing 300mbar negative pressure and capable of running continuously, (adequate cooling required) +5°C - +35°C. Existing pipework will not need replacing.

By reducing the negative pressure to 210mbar with suitable diameter pipework, this system is capable of running up to 6 surgeries.

Turbo-Smart 2V will operate from a 230V 13A supply.

Turbo Smart 2V delivers constant high airflow and vacuum, this machine automatically responds to match surgeries' demands.

Innovative, variable speed suction instantly adapts to any change in use.

Only one central separator is required rather than one per surgery, to minimise cleaning and maintenance. Cattani 'Turbo Smart 2V' is extremely reliable and easy to use, conveying performance data on a multi-function digital display.

Model / Code	Turbo Smart 2V 035135
Number of surgeries	2-6
Rated current	9 A
Rated kW	1.8 kW
Suction control	Variable speed suction
Maximum RPM	4800 RPM
Maximum negative pressure	300mbar
Sound level	69.8 dB(A)
Housing	NA
Operating environment	+ 5°C - + 35°C
Dimensions	Height: 650mm Depth: 360mm Width: 423mm Weight: 41 kgs
ISO-18 Amalgam retention system	98.1% retention rate at maximum flow of 18 l/min and minimum flow of 0.5 l/min
Standard complies to ISO 11143 - Dental Equipment Amalgam Separators	
Warranty	3 years