

MULTIPLE SURGERY ASPIRATORS

TURBO SMART Control The Power

Cattani 'Turbo Smart' is the powerhouse of suction systems, for simultaneous use by multiple surgeries. Delivering consistently high airflow and vacuum, this machine automatically responds to match surgeries' demands

Innovative, variable speed suction instantly adapts to any change in use. It speeds up when a suction tip is opened and when other surgeries are using the unit. When demand decreases, so does the machine's speed but without affecting performance.

Only one central separator is required rather than one per surgery, to minimise cleaning and maintenance. Cattani 'Turbo Smart' is extremely reliable and easy to use, conveying performance data on a multi-function digital display.

Model / Code	Model A 035115	Model B (inc 035500 software upgrade)
Number of surgeries	1-2	3-4 (parallel combined capability up to 16 surgeries)
Rated current	7 A	9 A
Rated kW	1.22 kW	1.5 kW
Suction control	Variable speed suction	
Maximum RPM	4800 RPM	6200 RPM
Maximum negative pressure	210mbar	
Sound level	68.2 dB(A)	69.6 dB(A)
Sound level with cover	67.3 dB(A)	68.7 dB(A)
Housing	Protected environment	
Operating environment	+ 5°C - + 35°C	
Dimensions	Height: 650mm Depth: 335mm Width: 400mm Weight: 36kgs	
ISO-18 Amalgam retention system	98.1% retention rate at maximum flow of 18 l/min and minimum flow of 0.5 l/min	
Standard complies to ISO 11143 - Dental Equipment Amalgam Separators		
Warranty	3 years	