NLX nano / NLX plus





NLX nano for high- and low-speed applications



NLX plus for high- and low-speed applications, including endodontics



Integrated power and control

Increase productivity with seamless integration

Equip your A-dec 500 and A-dec 300 with advanced high- and low-speed micromotor technology. With factory-integrated A-dec NLX nano and A-dec NLX plus, your A-dec delivery system gives you ready access for improved control, productivity and ease.

- Eliminate reaching, twisting and turning with an all-in-one delivery system
- Gain effortless control from your A-dec touchpad
- Program settings to match all your preferences
- · Monitor speed and switch between forward and reverse
- Simplify your environment, clear up space



Enhance the treatment room with a perfect fit

Place your micromotor instruments where they belong: on your delivery system. Co-developed for seamless operation, the combined NSK and A-dec solution enables you to focus on the patient instead of the equipment. The easy-touch control from the A-dec deluxe touchpad, the intuitive power display, and programmable memory for preferred settings increase productivity and promote better outcomes.

Experience a new generation of intuitive control

With applications ranging from high-speed gross reduction and low-speed precision finishing, to specialty functions such as endodontics, NLX nano and NLX plus help you achieve greater precision and treat patients more efficiently. Reduce fatigue: lightweight titanium, balanced handling ensure optimal ergonomics

Increase performance: full range of speed and all the power you need

Operate with precision and consistency: easy settings for activation, speed, torque, and direction

Do more in less time: activate the handpiece as you lift from the holder and view display as needed

Reduce cleanup and extend product life: autoclavable, maintenance-free, and brushless

Light up the procedure: LED daylight quality illumination gives you 32,000 lux for a clear view

Details	NLX nano	NLX plus
Applications	General restorative, prophylaxis	General restorative, prophylaxis, endodontics
Manufacturing Partner	NSK	NSK
A-dec Delivery System	A-dec 500/A-dec 300	A-dec 500/A-dec 300
Micromotor Controller*	NLX BA	NLX BA
Micromotor & Tubing	A-dec nano (factory) or NSK nano (in-field installation)	A-dec plus (factory) or NSK plus (in-field installation)
Contra-angle Compatibility	E-Style	E-Style
Motor Type	Sensorless DC Brushless Micromotor	Sensorless DC Brushless Micromotor
Speed	1,000 - 40,000 min ⁻¹	100 - 40,000 min ⁻¹
Torque	Max. 3.4 Ncm	Max. 4.0 Ncm, with Endo Function
Weight	61g (micromotor only)	70g (micromotor only)
Light	White LED 32,000 lux (micromotor only)	White LED 32,000 lux (micromotor only)

[*] An electric motor controller is a required component of the A-dec delivery system, which is configured at the time of your A-dec purchase.



A-dec Headquarters

2601 Crestview Drive Newberg, Oregon 97132 USA Tel: 1.800.547.1883 Within USA/CAN Tel: 1.503.538.7478 Outside USA/CAN a-dec.com.au Fax: 1.503.538.0276 a-dec.com

A-dec Australia

Tel: 1.800.225.010 Within AUS Tel: +61.(0)2.8332.4000 Outside AUS Tel: +44.(0).24.7635.0901 Outside UK

International Distribution Centers A-dec United Kingdom

Tel: 0800.ADEC.UK (2332.85) Within UK a-dec.co.uk

© 2012 A-dec Inc. All rights reserved. 85.0299.00/7-12/Rev A