

# CLEANI

Taking the work out of workflow

## PORTABLE AND COMPACT

Ideal chairside solution for any in-house workflows.

## BUILD PLATFORM WASH

Clean the build platform with prints all at once.

## DUAL TANKS

Significantly more washing space allows for heavy and fine rinses and prevents resin contamination.

## ACCELERATED WORKFLOW

90% effective cleaning process perfect for entry level to advanced users.



## PRODUCT FEATURES



### A Faster and Cleaner Workflow

CLEANI significantly reduces the time needed to clean prints while being able to clean up to 90% resin residue in one go. The new features and functions help dental clinics accelerate workflow by 40%.



### Alcohol Consumption Cut by Half

Alcohol lasts longer because the dual-tank design reduces cross-contamination and keeps the reservoirs cleaner. It is not necessary to repeatedly fill the tanks with alcohol after you print with a different resin.



### Wash Build Platform and Prints Together

CLEANI is built to accommodate both the prints and print platforms. The metal rods inside each tank can hold up most build and print platforms available on the market. This eliminates the need to first remove your prints from the print platform and then wash each piece separately.



### Vortex Minimizes Print Damage

The vortical motion creates gentle vibrations that minimize the chance of damage and breakage of prints during the wash cycle.



### Bigger Cleaning Space

Cleaning space is increased from the traditional one tank to two tanks, which allows for a heavy wash and fine wash for more prints.



### Safety Nets for Small Prints

No need to fish for prints anymore as the fine mesh nets in each tank catch small components and keep them from slipping through the cracks.



### Easy Maintenance and Cleaning

Removable tanks make maintenance and cleaning a much easier process. Instead of moving the entire machine, just simply remove the tanks and clean!



### Compact Design

A small washer with a big impact. CLEANI is ideal for dental clinics, as it's compact size allows you to take the entire cleaning process in-house.





## CLEANI

### Two tanks work better than one

CLEANI is designed with dual tanks for a two-stage cleaning process that allows for heavy and fine washes, with a workflow that's three times faster.

Ensure your 3D prints are cleaned properly before use and eliminate tedious, ineffective manual washes. Reduce and save on alcohol consumption with a stunning new workflow: CLEANI, the new dual-tank washer.

# ACKURETTA ECOSYSTEM

Envision a future where dental clinics have a solution to 3D print end-product restorations chairside with lower cost, and higher efficiency. From hardware, software, finishing solutions, materials, to training resources and documentation, Ackuretta keeps everything in-house and up-to-date.



## 3D Printing Materials

A complete and versatile resin portfolio designed specifically for dental applications. CURO resins produce high-strength, high-performance prints that are durable and offer enhanced patient comfort.



## Finishing Kit

Transform your 3D prints into safe, ready-to-use dental applications. The Ackuretta Finishing Kit effectively cleans and cures your prints to a professional grade finish.

## Software

A proprietary 3D slicing application, Alpha 3D is the most intuitive, UX driven, dental-focused, 3D slicing software on the market today.



**Alpha 3D**

## Other

All of this comes together to provide dental clinics and laboratories with a turnkey solution for digital dentistry



[sales@ackuretta.com](mailto:sales@ackuretta.com)



[ackuretta](https://www.ackuretta.com)



[ackuretta.com](https://www.ackuretta.com)