

SOL Breathe SOL into Your Creation


The perfect **SAME-DAY** chairside solution for **DIGITAL DENTISTRY.**

COMPACT, POWERFUL, AND ACCURATE, ready to print right out of the box.

Embrace SOL for your **DENTAL WORKFLOW** today, and start creating smiles for the future.







PRODUCT FEATURES



54 LEDs calibrated for extra curing consistency



150+ Prevalidated Materials



Selective curing Extended accuracy and lifespan



49µm pixel-size Super smooth surface finish



Instant access to resin updates



























APPLICATIONS



Dental Models
3 per print
50 min



Surgical Guides
10 per print
50 min



Temporary Crowns
55 per print
30 min



Denture Bases 4 per print 120 min



Denture Teeth
3 per print
45 min



Removable Partial Dentures
4 per print
60 min



Gingiva 45 per print 40 min



Impression Trays
2 per print
90 min



IBT 3 per print 80 min




Castable Crowns
35 per print
30 min



Mouth Guards
3 per print
70 min



Splints
3 per print
60 min



ACKURETTA ECOSYSTEM

Envision a future where dental clinics have a solution to 3D print end-product restorations chairside with lower cost, and higher efficiency. From hardware, software, finishing solutions, materials, to training resources and documentation, Ackuretta keeps everything in-house and up-to-date.



3D Printing Materials

A complete and versatile resin portfolio designed specifically for dental applications.



Finishing Kit

Transform your 3D prints into safe, ready-to-use dental applications.





