

Handpiece Special Offers

January - March 2023

TURBINES

T1 Boost Turbines

Perfection that fits in your hand



- Premium class LED turbines with ergonomic lightweight titanium sleeves
- Offering high power and optimum control at all times, the T1 range are extremely quiet during operation
- KaVo & Sirona fitting
- 3 year warranty



SIR-6689405
T1 Boost (KaVo Multiflex coupling)

RRP £1,514
Only £692

54% OFF

T2 Boost Turbines

Quality you can rely on



- You can rely on Dentsply Sirona T2 class LED turbines for any situation
- The continuous high performance, the ergonomic shape together with a surface finish for strong grip, these LED turbines made of stainless steel are guaranteed to satisfy
- KaVo, Sirona, W&H & NSK fitting
- 2 year warranty



SIR-6689637
T2 Boost (KaVo Multiflex coupling)

RRP £1,230
Only £527

57% OFF

T3 Boost Turbines

Powerful and economical



- Non optic
- Made of stainless steel, these non-fibre optic turbines have 26 watts of power and ceramic ball bearings to ensure efficient working in the long run
- Powerful lasting efficiency
- KaVo & Sirona fitting
- 1 year warranty



SIR-6689827
T3 Boost (KaVo Multiflex coupling)

RRP £903
Only £413

54% OFF

MASTERtorque Lux M9000L

KaVo's Premium Turbine



- Powerful - 23 watts
- Quiet - Only 57dB, vibration free
- Small head for improved access and patient comfort - Height 13.1mm
- Features Direct Stop Technology
- Four port spray
- Plasmatec coating for perfect grip
- Ceramic ball bearings
- Powerful illumination. Glass rod light conductor
- KaVo & Sirona fitting
- 2 year warranty



KAVO-M9000L
RRP £1,583
£950

Buy 3 get 1 FREE
(Only £713 each)

EXPERTtorque E680L



- Quiet, vibration-free operation (65 dB)
- 18 Watt power
- KaVo original ceramic ball bearings
- KaVo pushbutton chuck system (up to 30 N retention power)
- Compact head casing (height: 13.6 mm/Ø:12.5 mm)
- 4 port spray for optimal cooling
- Replaceable micro water filter
- KaVo MULTiflex coupling system
- Can be sterilised in an autoclave up to 135 °C
- KaVo & Sirona fitting
- 18 month warranty



KAVO-E680L
KAVO-E680L-DUOPK
RRP £1,188
£665

Or Buy 2 for £594 each

SMARTtorque S619L



- Quiet, vibration-free operation (65dB)
- 16 Watt power at only 2.8 bar
- KaVo pushbutton chuck system (30 N retention power)
- Compact head casing (height: 14.5 mm/Ø: 12.5 mm)
- 4 port spray for optimal cooling
- Can be sterilised in an autoclave up to 135 °C
- LUX glass rod fibre-optic
- KaVo fitting
- 12 month warranty




KAVO-S619L
RRP £727
£444

Buy 2 get 1 FREE
(Only £296 each)

Ti-Max Z900KL



- Standard Head
- Power 26 W
- Rotation 320,000-400,000 min-1
- Head Size Ø12.5 x H 13.1 mm
- Water Spray Quattro Spray
- Titanium Body Coating
- Cellular Glass Optics
- Ceramic Bearings
- Push Button Chuck
- Microfilter
- Fits Kavo, NSK, W&H, Bien Air, Sirona
- 3 years warranty




NSK-P1113001
RRP £1,068
£854

Buy 2 get 1 FREE
(Only £570 each)

S-Max M900L



- Standard Head
- Slim design with stainless steel body
- Speed: 325,000-430,000 min-1
- Head size: 12.1mm diameter x 13.3mm height
- 26W power
- Standard head
- NSK, Sirona, W&H, Bien-Air
- 2 years warranty




NSK-P1254001
RRP £747
£546

Buy 3 get 1 FREE
(Only £410 each)

S-Max M900KL



- Standard Head
- Slim design with stainless steel body
- Speed: 325,000-430,000 min-1
- Head size: 12.1mm diameter x 13.3mm height
- 26W power
- Standard head
- KaVo fitting only
- 2 years warranty



NSK-P1258001
RRP £747
£393

DMI exclusive 47% OFF

Pana-Max 2



- Non optic
- Rotation 350,000-450,000 min-1
- Head Size Ø11.2 x H 13.4 mm
- Water Spray Single Spray
- Stainless Steel
- Ceramic Bearings
- Push Button Chuck
- Interchangeable Cartridge
- Fits KaVo & NSK
- 1 year warranty



NSK-P1212001
RRP £464

Buy 1 get 1 FREE
(Only £232 each)

Synea Vision TK-98L

- Ring LED+,
- 5-hole spray,
- Head size Ø 11.5 mm,
- 21 W, 360,000 rpm,
- FG burs Ø 1.6 mm for bur lengths up to 25mm
- Ceramic ball bearings,
- Coupling system Roto Quick



RRP £1,563
£859

WH-30019000

Synea Fusion TG-98LM

- 4 port spray
- Head size Ø 11.5mm,
- 19 W, 360,000 rpm,
- FG burs Ø 1.6mm for bur lengths up to 25mm
- Compact glass rod, ceramic ball bearings,
- Coupling system Multiflex, NSK & Rotoquick



RRP £1,069
£588

WH-30003000

Offers cannot be combined with any other promotions. Excludes VAT E&OE.

Call: 028 9260 1000

Web: www.dmi.co.uk

Email: info@dmico.uk



CONTRA ANGLES 1:1 RATIO (BLUE)

T1 C40L



- Longer service life thanks to Diamond like Carbon coated ball bearings
- Optimal illumination for an excellent view
- 1:1 ratio has single spray
- 3 year warranty

SIR-3328437 Model: C40L
1:1 Ratio (Blue)

RRP £937

Now £499

47% OFF



T2 A40L



- Robust ceramic ball bearings
- Pleasant and easy to grip thanks to its titanium sleeves
- Optimal illumination for great visibility during treatments
- 1:1 ratio has single spray
- 2 year warranty

SIR-6376185 Model: A40L
1:1 Ratio (Blue)

RRP £783

Now £469

40% OFF



T3 E40



- The T3 is a sturdy, economical contra angle which will give many years of high performance
- Sturdy, titanium sleeve
- Ceramic ball bearings
- Non optic
- Internal water
- 1 year warranty

SIR-6392687 Model: E40
1:1 Ratio (Blue)

RRP £572

Now £383

33% OFF



MASTERmatic M20L & L68



- Transmission 1:1
- Speed range up to 40,000 rpm
- Optic
- Internal air/water coolant
- Works on electric or air motors that have a universal "E" type connection
- Plasmatic Anti-Slip Finish
- Improved angulation
- Flexibility of use when paired with different INTRA LUX Heads
- 2 year warranty

KAVO-M20L

RRP £1,201

£758

Buy 3 get 1 FREE
(Only £569 each)



EXPERTmatic E20L



- Quiet
- Transmission 1:1
- Cellular optic
- For latch type burs
- Speed range up to 40,000 rpm (electric motor) and 20,000 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip™
- Universal "E" type connection
- 18 month warranty

KAVO-E20L-DUOPK
KAVO-E20L

RRP £958

£587

Or Buy 2 for £479 each



SMARTmatic S20

- Non optic, no water spray
- Ceramic sliders within the clamping system for enduring and secure tool holding
- Rust-free stainless steel for long-lasting durability
- Lightweight for tireless working
- Optimum balance for precise treatment
- 1:1 Ratio
- 1 year warranty

KAVO-S20

RRP £236

£170



Ti-Max Z25L



- Anti Heat System
- Cellular Glass Optics
- Clean Head System
- Push Button Chuck
- Microfilter
- Ratio: 1:1
- Head Size: H13.7 x Ø8.9
- Single Spray
- Titanium
- Max Speed: 40,000 min-1
- 3 year warranty

NSK-C1038001

RRP £1,128

£903

Buy 2 get 1 FREE
(Only £602 each)



S-Max M25L



- Ratio: 1:1
- Optic
- For CA burs (Ø2.35)
- Water Spray: Single Spray
- Body Material: Stainless Steel
- Max Speed: 40,000 min-1
- 2 years warranty

NSK-C1024001

RRP £745

£596

Buy 3 get 1 FREE
(Only £447 each)



FX25M



- Non optic. No water
- Rotation 350,000-450,000 min-1
- Head Size Ø11.2 x H 13.4 mm
- Stainless Steel
- Ceramic Bearings
- Push Button Chuck
- Interchangeable Cartridge
- 1 year warranty

NSK-C1052001

£251

Buy 1 get 1 FREE
(Only £126 each)



Synea Vision WG-56LT



- Contra-angle handpiece 1:1
- Single Spray
- Optic
- Head size Ø 9.5 mm
- Push-button chuck for contra-angle bur with Ø 2.35 mm
- Compact glass rod
- 2 year warranty

WH-30009000

RRP £1,069
£588



Alegria WE-56



- Non optic/no spray
- Contra-angle handpiece 1:1
- Push-button chuck system
- For contra-angle burs with Ø 2.35 mm (ISO 1797-1)
- Drive speed: 40,000 min-1
- Outer diameter of the sheath: 20 mm
- Head size: Ø 10.1 mm
- Head height (with 19 mm bur): 20.8 mm

WH-30116000

RRP £312
£120



SPEED INCREASING CONTRA ANGLES 1:5 (RED)

T1 C200L



- Longer service life thanks to Diamond like Carbon coated ball bearings
- Optimal illumination for an excellent view
- 1:5 ratio has four hole spray
- 3 year warranty



SIR-3328452
Model: C200L 1:5 Ratio (Red)

RRP £1,505
Now **£804**

47% OFF

T2 A200L



- Robust ceramic ball bearings
- Pleasant and easy to grip thanks to its titanium sleeves
- Optimal illumination for great visibility during treatments
- 1:5 ratio has four hole spray
- 2 year warranty



SIR-6376177
Model: A200L 1:5 Ratio (Red)

RRP £1,058
Now **£636**

40% OFF

T3 E200



- Contra-Angle Red
- Ratio 1:5
- Non optic
- Spray 4N
- Ceramic bearings
- Head Diameter 8.7
- Head Height 14.6
- RPM 200,000
- 1 year warranty



SIR-6392679

RRP £847
Now **£563**

34% OFF

MASTERmatic M25L

- Speed increasing handpiece ideal for crown and veneer preparation
- 1:5 speed increasing ratio
- Extremely quiet - Only 55dB
- Features Direct Stop Technology
- No suck back. Replaceable spray microfilter
- Three port spray
- Plasmatec coating for perfect grip
- Ceramic ball bearings
- Powerful illumination
- 2 year warranty



KAVO-M25L
RRP £1,789
£1,026

Buy 3 get 1 FREE
(Only £769 each)

EXPERTmatic E25L

- Speed-increasing handpiece transmission 1:5 (red)
- Fast and precise working thanks to even smaller head and ceramic bearings
- Best reliability thanks to internally conducted spray water and air, 3-port spray, and replaceable water filter
- Perfect lighting thanks to the 25,000 Lux glass rod light conductor
- Better Grip
- 18 month warranty



KAVO-E25L-DUOPK
KAVO-E25L
RRP £1,314
£788

Buy 2 for £657 each

Ti-Max Z95L Switch

- Optic contra-angle
- 1:5 Increasing
- Quattro spray
- Ceramic ball-bearings
- World first switch between water and spray on handpiece
- 3 year Warranty



NSK-C1034001SW
£1,112

Buy 2 get 1 FREE
(Only £741 each)

S-Max M95L

- 1:5 Increasing
- Max Speed 200,000 min-1
- Head Size Ø9.9 x H 14.5 mm
- For FG Burs (Ø1.6)
- Body Material Stainless Steel
- Water Spray Quattro Spray
- Optic
- 2 year warranty



NSK-C1023001
£735

Buy 3 get 1 FREE
(Only £552 each)

S-Max M95

- Non Optic
- 1:5 Increasing
- Max Speed 200,000 min-1
- Head Size Ø9.9 x H 14.5 mm
- For FG Burs (Ø1.6)
- Body Material Stainless Steel
- Water Spray Quattro Spray
- 2 year warranty



NSK-C1026001
£596

Buy 3 get 1 FREE
(Only £447 each)

EVO.15 1:5 L



- Standard or Micro Series
- Small head (10% smaller)
- CoolTouch+™ heat-arresting technology
- Shockproof head
- 55 dBA
- Multi-strand optical glass technology
- Up to 3 years warranty

RRP £1,227
£773



BIEN-1600940-001 Micro Series
BIEN-1600941-001 Standard

CA 1:5 L



- Standard or Micro Series
- Multi-strand optical glass technology
- 58 dBA
- Cool Touch™
- Up to 3 Years Warranty

RRP £1,023
£550



BIEN-1600386-001 Standard
BIEN-1600690-001 Micro Series

Or Buy 2 for £1,054

AIR MOTORS

AM1006 with LED

- Chrome coating
- LED Light
- 5,500 K colour temperature
- 20,000 rpm non-load speed
- Left and right direction of rotation
- Internal spray
- 6-pin tubing
- 1 year warranty

MKDE-AM1006

£339

MK-dent
GERMANY



181L with LED

- Speed range up to 20,000 rpm
- Forward / reverse speed adjustment
- Internal air / water coolant
- 6-pin connection
- Fits all attachments optic or non-optic with Universal "E" type connections
- 2 year warranty

KAVO-181L

£961

KAVO
Dental Excellence



AM-25RM

- Non optic
- No spray
- New integrated 360° rotation
- Mid west connection
- 1 year warranty

W&H-30108000

£218

W&H



AM0004

- Max. 20 000rpm
- No spray
- Standard INTRAmatic® connection
- Suitable for all common instruments
- Non fibre optic
- 1 year warranty

MKDE-AM0004

£175

MK-dent
GERMANY



STRAIGHT HANDPIECES 1:1 RATIO

Allegra HE-43

- For handpiece and contra-angle burs with Ø 2.35 mm (ISO 1797-1)
- Drive speed: 40,000 min-1
- Outer diameter of the sheath: 20 mm
- Non optic
- No spray
- 1 year warranty

WH-30119000

£115

W&H



T4 Model BH40

- Stainless steel
- Bur fastening via push button
- For use with burs with a 2.35mm diameter stem
- Max. speed: 40,000rpm
- No spray
- No optic
- 1-year warranty

SIR-6401843

£151

Dentsply
Sirona



Smartmatic S10

- 1:1 (blue handpiece)
- Main application: Restorative
- Max. drive operating speed: 40,000
- Non optic
- No spray
- 1 year warranty

Kavo-S10

£163

KAVO
Dental Excellence



FX65M

NSK

- 1:1 Straight
- No spray
- Max Speed 40,000 min-1
- No spray
- Non optic
- 1 year warranty

NSK-H1014001

£251



Buy 1 get
1 FREE
(Only £126 each)

COUPLINGS FOR KAVO TYPE TURBINES

KaVo Multiflex 465LED

- LED Optic with water control
- With 4-hole standard connections
- Rotary ring spray regulation with multiple indent locking positions
- Non-return valve (water)
- 1 year warranty

Kavo-465LED

£317

KAVO
Dental Excellence



KaVo Multiflex 465LRN

- Optic with water control
- With 4-hole standard connections
- New lamp holder
- Rotary ring spray regulation with multiple indentlocking positions
- Non-return valve (water)
- 1 year warranty

Kavo-465LRN

£269

KAVO
Dental Excellence



KaVo Multiflex 465RN

- With 4-hole standard connections
- Rotary ring spray regulation with multiple indent locking positions
- Non-return valve (water)
- 1 year warranty

Kavo-465RN

£173

KAVO
Dental Excellence



KaVo Multiflex 460LED

- LED Optic
- No water control
- 1 year warranty

Kavo-460LED

£175

KAVO
Dental Excellence



KaVo Multiflex 460E

- Non optic
- No water control
- 1 year warranty

Kavo-460E

£118

KAVO
Dental Excellence



MK DENT Multiflex

- 6 hole quick connector with LED
- With water control
- 1 year warranty

MK-dent
GERMANY

MKDE-QC6016KW

£179



MK DENT Multiflex

- 6 hole quick connector with LED
- No water control
- 1 year warranty

MKDE-QC6016K

£149

MK-dent
GERMANY



COUPLINGS FOR W&H TYPE TURBINES

RQ-34

W&H

- Optic with water control
- Flexible
- Free 360° instrument rotation
- Continuous regulation of the spray water supply
- Roto Quick non-retraction valve
- 1 year warranty

WH-10403400

£168



RQ-24

W&H

- Optic
- No water control
- Roto Quick non-retraction valve
- 1 year warranty

WH-10402400

£152



RQ-14

W&H

- Non optic
- Water control
- Free 360° instrument rotation
- Continuous regulation of the spray water supply
- Roto Quick non-retraction valve
- 1 year warranty

WH-10401400

£128



RQ-04

W&H

- Non Optic
- No water control
- 1 year warranty

WH-10400400

£96



Call: 028 9260 1000

Web: www.dmi.co.uk

Email: info@dmico.uk

HANDPIECE MAINTENANCE/AIR POLISHERS

Sirona LED coupling

- Optic
- 1 year warranty

SIR-6319706



Dentsply Sirona

£324

Coupling for NSK - PTL CL

- LED optic coupling
- No water control
- 1 year warranty

NSK-P1001600



NSK

£175

QUATTROcare Plus 2124A

The cost-efficient solution for optimum care and long-term maintenance of handpieces.

- Cleans 3 slow and 1 fast handpiece or 2 fast and 2 slow handpieces
- Only one-minute care cycle
- Program-controlled, automated care process
- Optimum lubrication, thanks to automatic, outlet-specific, metered dosing
- Easy handling
- No leakage of excess residual oil from handpiece
- Modern and ergonomic design

KAVO-1-008-3805

KAVO
Dental Excellence



RRP £2,340 **£1,638**

PROPHYflex 4 Air Polisher Box

- Easy to use Two handpieces in one – perfect for supragingival and subgingival applications - Whether for prophylaxis or periodontology, fissure sealing, aesthetic dentistry, prosthetics or orthodontics, it can be used on any teeth, any patient and almost any indication
- PROPHYflex 4 is available in three fresh and exciting colors, each available as a MULTiflex coupling or Sirona coupling option

KAVO
Dental Excellence

Kavo-1-013-5078
Blue
Kavo-1-013-5079
Lime
Kavo-1-013-5080
Flamingo



Each PROPHYflex Box includes:

- 1x PROPHYflex 4 Air Polisher
- 1x PROPHYpearls neutral (80 sticks)
- 1x PROPHYpearls selection (80 sticks)
- 1x Kerr OptiView Standard Kit (2 pcs)

£1,128

Perio Box

Enables sub gingival cleaning
Contains: 1 perio kit, 1 perio powder (4 x 100g),
1 Kerr perio probe (2pcs) Kavo-1-013-5081

£280

RONDOflex plus 360

- Highly efficient, clean work, light-weight and accurate
- Gentle treatment without drilling, without injection, without noise, without vibrations. Pleasant for your patient, particularly for sensitive patients and children
- Incomparable flexibility through problem-free use at various workplaces thanks to MULTiflex coupling: simply push on and start.
- Includes 2 cannulae.

RRP £1,847

KAVO-1-002-2179



KAVO
Dental Excellence

**Now Only
£1,319**

NSK iCare

- Cleans 4 Handpieces
- Internal cleaning and lubrication
- Chuck cleaning function ensures that the chuck is thoroughly cleaned and lubricated
- Purging select mode ensures complete removal of excess oil
- Integrated large oil tank (1.2 litre) reduces burden of refill
- 28cm (W) X 24cm (D) X 35cm (H)
- Air Line Required

NSK-Y1002787

NSK



£1,286

Prophy-Mate Neo

- Air Powered Tooth Polishing System
- Twin nozzle
- Comfort design
- Twin joints rotate smoothly 360° degrees for easy operability
- One touch operation
- Selectable 60° and 80° nozzles
- The handpiece can be washed in the thermoisinfectior
- Autoclavable up to 135°C
- Fits KaVo, NSK, Sirona, W&H, Bien Air

NSK



£990

NSK-Y135030

Universal Oil Spray

- High-quality universal care for all KaVo contra-angles, handpieces, turbines, heads and air motors
- Ensures a long life for all KaVo instruments
- Environmentally friendly propellant
- 500ml

KAVO
Dental Excellence



KAVO-0-411-9640

£27.49

Quattrocare Plus Oil Spray

- Designed to clean, lubricate and purge your instruments
- Maintenance spray for KaVo turbines, air motors, straight and contra-angle handpieces
- Only for use in the QUATTROcare handpiece oiler
- 500ml

KAVO
Dental Excellence



KAVO-1-011-5720

£27.49

Oil Nozzle Intra (Slow handpieces)

- Handpiece lubricating oil spray nozzle for slow handpieces

KAVO
Dental Excellence



KAVO-0-411-9911

£8.56

Oil Nozzle Multiflex (Fast handpieces)

- Handpiece lubricating oil spray nozzle for fast handpieces

KAVO
Dental Excellence



KAVO-0-411-9921

£8.56

Oil Spray

- For use with all types of handpiece
- Non-toxic and non-allergenic
- No odour or taste
- Good wetting capacity to cover all internal moving parts
- Heat resistant and does not evaporate
- 500ml

PEGA-100111

£11.10



T1 Oil Spray

- Universal cleaning and lubricating oil for the preservation and disinfection of dental instruments
- 250ml

Dentsply Sirona



SIR-5901665

£39

F1 Service Oil Spray MD

- Superior lubrication
- Fully sterilizable, absolutely resistant to hydrolysis
- F1 Oil works with turbines, straight and contra-angle handpieces, surgical instruments, air scalars and air motors
- 400ml

WH-10940021

£36.95



Oil Spray Nozzle - Fixed Connection

- Spray Cap With Nozzle
- Chuck system for turbines and contra-angle handpieces as well as turbine handpieces and air motors with fixed connection

WH-02036100

£5.95



Oil Spray Nozzle - Slow Handpieces

- Spray cap for W&H F-1 service oil
- For slow handpieces

WH-02038200

£5.95



Oil Spray Nozzle - Fast Rotoquick Connection

- Spray cap with nozzle for all W&H turbines, air motors and scalars with ROTO QUICK connection

WH-02229200

£15.92



Pana Spray Oil

- PANA SPRAY Plus is a specially formulated, highly effective general handpiece and air motor lubricant.

NSK



NSK-Z182100

£29.43

iCare Maintenance Oil

- Specially formulated maintenance oil provides for the highest levels of cleaning & lubricating of handpieces
- 1000ml

NSK



NSK-Z016130

£84.39

Woodpecker O - Star

The Ultimate Curing Expert

- Performs light curing in 1 second
- Longer wavelength: 385 - 515 nm
- It includes a light intensity meter which is located in the charging base
- It has 7 different work modes: normal, high intensity, turbo, ortho, soft, pulse and verification
- Light intensity: up to 3000mw/cm²
- Small, rotating head allows you to reach all angles of the mouth and reach the back teeth without difficulty

WOOD-OSTAR

£495

**NEW
Premium Light**

WOODPECKER



Woodpecker iLED

Curing Light

- Ergonomic curing light with 360° rotating head
- Fast - Cures up to 2mm in just one second
- Light intensity up to 2300mw/cm² with constant power output
- Power input: AC100V-240V
- Light output intensity: 1000 mW/cm²~2500mW/cm²
- Time Sets: 5-20 Seconds
- Operation Mode: turbo mode, normal mode

£250

WOOD-ILED

WOODPECKER



Woodpecker LED.H ORTHO

Specially designed for Orthodontics

- Specially designed for orthodontics. In just 3 seconds the brace solidifies
- Short, powerful cycles without overheating
- Constant work without losing intensity, for a complete curing of the resins
- Light and easy to use, with a rechargeable battery
- The new bevelled light guide design can easily reach any tooth and radiate vertically at any angle
- Includes a light intensity meter

WOOD-LED.H.ORTHO

£195

WOODPECKER



Woodpecker LED.H

Specially designed for Orthodontics

- Woodpecker LED-H Light Curing Lamp provides reliable light curing of any type of material
- It is a light weight lamp, which will allow you to use it during a large part of your professional daily activity, without causing fatigue
- Three working modes: Normal, ramp or intermittent
- High-capacity battery with charge duration for 200 continuous light-curing cycles

WOOD-LED.H

£160

WOODPECKER



Sine Wave Reservoir Scaler

- Ultrasonic cleaning device with a handpiece type EMS without light
- Dismountable and sterilizable table-top model
- Modern design with touch control and waterproof
- Scaler, perio and endo function
- Modulates sinus waves, allowing the best treatment with maximum comfort

Includes:

- 7 tips: 1 x G4 supragingival scaler tip
- 1 x G12 supragingival scaler tip
- 1 x P1 Periodontal scaler tip
- 1 x P3 Periodontal scaler tip
- 1 x P8 Periodontal scaler tip
- 2 x E6 Endodontic scaler tip

WOOD-U600

£595

WOODPECKER



Ultrasonic UDS-N2 Scaler

- Built in ultrasonic scaler (compatible with A-dec and Belmont units)
- The handpiece is detachable and can be autoclaved to 135 °
- The UDS N2 instrument is compatible with the EMS system compatible tip and inserts

WOODPECKER



£395

Price includes installation

WOOD-UDS-N2

DTE-D7 Piezo Scaler

Satelec Type

- Assembled with excellent piezo-ceramic, the handpiece has high and lasting power and generates little heat.
- Eliminates dental calculus, bacterial plaque and tea scales on the teeth quickly without hurting the gum or the enamel.
- Light weight, easy to operate.
- Compatible with Satelec hand piece and tips.



£425

Replacement Handpieces

- Woodpecker/EMS Type WOOD-HW-3H £145
- Woodpecker U600 Handpiece WOOD-HW-6 £145
- Woodpecker/EMS Type LED WOOD-HW-5L £160
- DTE/Acteon Satelec Type LED DTE-HD-7H £145

Replacement Wrenches

- Woodpecker/EMS Type WOOD-TW-IL £10.50
- Sirona Type WOOD-TS-IL £13
- DTE/Acteon Satelec Type DTE-TD-IL £10.50
- Satelec Type DTE-TD-2L £13

Replacement Tips

- Woodpecker/EMS
- Perio P1 WOOD-P1
- Perio P3 WOOD-P3
- Perio P4 WOOD-P4
- Scaling G1 WOOD-G1

- DTE/Acteon/Satelec
- Perio PD1 DTE-PD1
- Perio PD3 DTE-PD3
- Scaling GD1 DTE-GD1
- Scaling GD4 DTE-GD4

£18